

Monitor your Miri® incubator

Connect up to eight (8) Miri incubators into a PC or a tablet PC to monitor the performance of your incubators. The logging software will document parameters and optionally send e-mail alarms.

You can order from us a set of tablet PC and a holder for data-logging of your Miri® or you can buy your own tablet with the following specifications:

Minimum system requirements for Miri with SafeSens

- Intel Core 2 Duo or AMD Athlon X2 at 2.4 GHz processor
- 4GB RAM
- 15 Gb Hard Disk space
- Integrated video card
- Monitor with resolution of 1024 x 768
- Windows 7 Pro/ 8 Pro/ 10 Pro operating system with 64 Bit Architecture
- USB 3.0 port for each connected device

Minimum system requirement for Standard Miri

- Windows 7 or Higher
- USB 2.0 port for each connected device
- Tablet PC or desktop PC



ORDER NOW!

Item Code	Description
I320140	PC Holder for Miri®
I320141	Tablet PC for Data Logging of MRI-6A10
I320142	Set of Tablet PC and Holder For Data Logging of MRI-6A10

Esco Micro Pte. Ltd.

21 Changi South Street I
Singapore 486 777

Tel +65 6542 0833
Fax +65 6542 6920

medical@escoglobal.com
www.medical.escoglobal.com

Esco Global Offices: Bahrain | Bangladesh | China | Denmark | Germany | India | Indonesia | Italy | Japan | Malaysia | Philippines | Russia | Singapore South Africa | South Korea | Thailand | United Kingdom | USA | Vietnam

