

# SAFE Sens Continuous pH Monitoring

Worry-free pH monitoring of your Esco incubators

SAFE Sens (Sterile, Automated Fluorescopic Evaluation) offers fast, effective, and non-invasive continuous pH monitoring product for *in vitro* fertilization (IVF) procedures.

The SAFE Sens technology employs an optical fluorescent measurement technology, used in combination with disposable sensors, which accurately and reliably monitors the pH of small volumes of fluids such as the media used in IVF.

Order SAFE Sens with your new **Miri®** and **Mini Miri®\***



\*Factory-installed. No option yet to install on-site for previously ordered Miri.



## Continuous pH measurement

- Reading and recording every 30 minutes (default setting - adjustable).
- Single use sensor probe for up to seven (7) days of pH readings.



## Data-Logging System

- Data Logging and user alarms.
- Each TrakStation® can be connected to multiple incubators.



## Easy to implement

- Easy to align (no buffers, no hassles).
- Easy to use and maintain.



## Compact and Efficient

- No more unnecessary openings of your incubator for spot pH measurement.
- Only requires 100 µL of media + 50 µL of oil.